

BRENMOOR

brenmoor.com



Economical single-use tourniquet **Tournistrip[®]eco**



Tournistrip[®]eco

Comfortable, easy to apply,
single-use tourniquet.

Designed by clinicians for ease of use and a better experience for patients, **Tournistrip[®]eco** allows effortless conversion from reusable devices to single-use, avoiding a potential transfer vehicle for the spread of nosocomial infections, including MRSA.

Tournistrip[®]eco looks and operates like a conventional tourniquet. Its unique design and construction means the tourniquet is truly single use. It is easy to tension, adjust and can be repositioned on the same patient up to 5 times. Release and removal is single-handed.

- Design - aids infection prevention and prevents cross infection
- Fitting - release with a single hand, minimising contact
- Adjustment - easy to tension, adjust and reposition
- Flexible - join multiple tourniquets for larger/bariatric patients
- Comfortable - doesn't pinch skin like 'rubber band' designs
- Sensitive - 100% latex free
- Proven - clinically trialled with 226 patients and 20 phlebotomists



ADJUSTABLE
LENGTH

SELF ADHESIVE
FASTENER

LATEX FREE
CONSTRUCTION

COMFORTABLE
FOR PATIENT

EASY TO FIT
TO PATIENT

OTHER COLOURS
AVAILABLE

SINGLE-USE

For more information visit brenmoor.com or email info@brenmoor.com

BRENMOOR

brenmoor.com

Economical single-use tourniquet Tournistrip[®]eco



SINGLE USE

Unique design and construction to aid infection prevention and prevent cross infection.

EASY TO USE

Easy to tension and adjust. Can be repositioned on the same patient up to 5 times.

SINGLE-HANDED RELEASE

Can be released with a single hand, minimising contact.

LATEX FREE

Tournistrip[®]eco is 100% Latex free.

COMFORT

Tournistrip[®]eco doesn't pinch the skin like "rubber band" designs.

FLEXIBLE

Two Tournistrip[®]eco can be joined together for larger arms/bariatric patients.

CLINICALLY PROVEN

Tournistrip[®]eco is the only medical tourniquet to have been clinically trialled, involving two London Teaching Hospitals with 226 patients and 20 phlebotomists.

Tournistrip dispenser with refill

Dispenser can be wall mounted, adhered to a flat surface or fixed to a medical grade trolley. When refilling the dispenser is unlocked via the front latch. The refill pack can then be inserted and the tourniquet pushed through the aperture.

ORDER CODE BM1081 ●
Orange

Refill pack of 150 Tournistrip[®]eco

Remove from liner and place the refill inside the dispenser on the centrally located holder. Unwind the first tourniquet on the strip and place through the aperture of the Tournistrip dispenser.


SPECIFICATIONS				
Box qty 900	Roll detail 6 x rolls of 150	Size 510 x 30mm	Core size 50mm	Outer diameter 120mm

ORDER CODE BM1078 ○
White

DIMENSIONS >>>

BRENMOOR

DESIGNED & MANUFACTURED IN THE UK



For more information visit brenmoor.com or email info@brenmoor.com